

Gladfelter Lab Protocols
Immunofluorescent Staining in *Ashbya*

- Grow a 50 mL culture of *Ashbya* at 30°C for 15 hours. Add 5.5 mL 3.7% formaldehyde and shake at 30°C for 1 hour.
- Prepare polylysine slides: Add 10 μ L polylysine to each well for 5 minutes (aliquots stored in common stock -20°C), aspirate and let air dry, wash with 10 μ L water, aspirate and let air dry.
- Divide culture into 15 mL conicals (cells do not pellet well in 50 mL conicals) and spin in benchtop centrifuge at 3000 rpm for 5 minutes.
- Remove supernatant (discard in formaldehyde waste in fume hood). Pool cells into one 15 mL conical and add 10 mL PBS. Spin at 3000 rpm for 5 minutes and remove supernatant.
- Add 1 mL Solution A and transfer cells to a 1.5 mL Eppendorf tube.
 - If you have many cells, you can add 2 mL Solution A and reserve 1 mL of cells at 4°C for up to a week.
- Add 100 μ L of 15 mg/mL Zymolyase (aliquots stored in common stock -20°C). Incubate at 37°C, gently inverting every 5 minutes. Check cells under microscope every 10 minutes until 70-80% of cells are phase dark. Wild type typically takes 30-40 minutes, but strains vary considerably.
 - Note: Zymolyase is a fragile enzyme, when thawing do not vigorously shake or vortex the tube.
- Centrifuge at 2600 rpm for 2 minutes, remove supernatant. Wash 2 times with 1 mL solution A and resuspend in 500 μ L Solution A.
 - Digested cells are fragile, take care while washing or they explode.
 - Final resuspension volume can be adjusted depending on cell density.
- Apply 20 μ L cells onto polylysine treated well. Allow cells to settle for 15 minutes and aspirate. Let air dry completely.
- Wash cells 2X with 20 μ L PBS.
- Block with 10 μ L PBS + 1 mg/mL BSA (IgG free, stored at 4°C) for 30 minutes.
- Wash 2X with 20 μ L PBS.
- Add 10 μ L primary antibody diluted in PBS+BSA and place slide in humid chamber. Incubate at 4°C overnight. (antibodies stored in common stock -20°C)
 - Rat α -alpha tubulin use 1:50
 - Rabbit α -Cdc11 use 1:100

- Wash 5X with 20 μ L PBS+BSA
- Incubate 1 hour in 1:200 secondary antibody + Hoechst (1:500) in PBS+BSA
 - o AlexaFluor 488 or 568 goat α -rat and AlexaFluor 568 rat α -rabbit use 1:200
- Wash 5X with 20 μ L PBS+BSA
- Add 7 μ L Prolong Gold mounting medium (aliquots stored in common stock -20°C).
Apply coverslip, let settle for 30 minutes with a kimwipe on top and weight applied, and seal edges with nail polish.

Solution A (100 mL)

10 mL 1 M Potassium phosphate buffer pH 7.5

60 mL 2 M Sorbitol

30 mL H₂O

filter sterilize