

PCR: NAT1 deletion cassette (ClonNAT)

Mix

Stock for 6 reactions à 100 µl (for 1 deletion):

| | |
|--------------------------|--|
| 209 µl | H ₂ O |
| 120 µl | 5x PCR Buffer (Roche: GC rich kit) |
| 60 µl | Resolution solution (Roche: GC rich kit) |
| 60 µl | dNTPs (Stock: 2 mM) |
| 25 µl | MgCl ₂ |
| 6 µl (more if low conc.) | pUC19NATPS template (100 – 300 ng) |
| 30 µl | Primer NS1 (Stock: 10 µM = 10 pmol/µl) |
| 30 µl | Primer NS2 |

Aliquot 90 µl to 6 PCR tubes, heat to 95 °C and add 10 µl polymerase mix
Polymerase Mix: **3.5 µl Poly** (Roche: GC rich kit) + **66.5 µl H₂O**

Cycle Parameters

| | | |
|----|--------------|--------------------------|
| 1. | 95 °C | 5' 00" |
| | | Add 10 µl polymerase mix |
| 2. | 95 °C | 0' 30" |
| 3. | 53 °C | 0' 45" |
| 4. | 72 °C | 2' 35" |
| | GOTO 2. | REP 30 x |
| 5. | 72 °C | 30' 00" |
| 6. | 4 °C | ∞ |

Collect the 6 reactions, precipitate, redissolve in 75 µl H₂O and transform 200 µl mycelium.

Complications: The buffers of the GC rich kit lead to high salt in the DNA pellet, which would reduce efficiency of electroporation.

⇒ use 50 % EtOH for the washing steps of DNA precipitation.

⇒ redissolve in 75 µl instead of 50 µl H₂O